

Product Datasheet

RNF123 Antibody (PerCP/Cy7) (orb1816250)

Description	RNF123 Antibody (PerCP/Cy7)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PerCP/Cy7
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human RNF123 (1-100/1314aa)
Target	RNF123
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5XPI4
Expiration Date	12 months from date of receipt.