

Product Datasheet

APP Antibody (Biotin) [Out of stock] (orb1816257)

Description	APP Antibody (Biotin) [Out of stock]
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human beta-Amyloid (1-40) (1-40/42aa)
Target	APP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P05067
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.