

Product Datasheet

CALCA Antibody (Biotin) (orb1816357)

Description	CALCA Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Calcitonin Protein (26-141/141aa)
Target	CALCA
Preservatives	0.01M PBS(pH7.4)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01258
Expiration Date	12 months from date of receipt.