

Product Datasheet

169565-17-7 Antibody (Biotin) (orb1816532)

Description	169565-17-7 Antibody (Biotin)
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Target	169565-17-7
Preservatives	0.01M TBS(pH7.4)
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
CAS Number	169565-17-7
Expiration Date	12 months from date of receipt.