

Product Datasheet

ALPP Antibody (orb1816704)

Description	ALPP Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Placental alkaline phosphatase
Target	ALPP
Preservatives	Liquid in PBS containing 50% Glycerol, 0.5% rAlbumin and 0.02% Proclin300.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1/kappa
Clonality	Monoclonal
MW	53 kDa
Uniprot ID	P05187
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.