

## Product Datasheet

### RPA2 Antibody (orb1816726)

|                            |                                                                                    |
|----------------------------|------------------------------------------------------------------------------------|
| <b>Description</b>         | RPA2 Antibody                                                                      |
| <b>Species/Host</b>        | Mouse                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                       |
| <b>Tested Applications</b> | IHC-P, WB                                                                          |
| <b>Immunogen</b>           | Recombinant human RPA2 Protein                                                     |
| <b>Target</b>              | RPA2                                                                               |
| <b>Preservatives</b>       | PBS with 0.02% Proclin300                                                          |
| <b>Concentration</b>       | 1mg/ml                                                                             |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.<br>Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                              |
| <b>Isotype</b>             | IgG1,K                                                                             |
| <b>Clonality</b>           | Monoclonal                                                                         |
| <b>MW</b>                  | 29 kDa                                                                             |
| <b>Uniprot ID</b>          | <b>P15927</b>                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                    |