

Product Datasheet

NES Antibody (orb1816769)

Description	NES Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Nestin
Target	NES
Preservatives	PBS, Glycerol, rAlbumin
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a κ
Clonality	Monoclonal
MW	178 kDa
Uniprot ID	P48681
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.