

Product Datasheet

GZMB Antibody (orb1816787)

Description	GZMB Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human GZMB
Target	GZMB
Preservatives	PBS, Glycerol, rAlbumin
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a, κ
Clonality	Recombinant
MW	25 kDa
Uniprot ID	P10144
Expiration Date	12 months from date of receipt.