

Product Datasheet

LAMP1 Antibody (orb1816814)

Description	LAMP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant human LAMP1 protein (29-382/417aa)
Target	LAMP1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11279
Dilution Range	IHC-P=1:100-500, WB=1:500-1:2000, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.