

Product Datasheet

RHOJ Antibody (orb1816849)

Description	RHOJ Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human RHOJ (61-160/214aa)
Target	RHOJ
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	Q9H4E5
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.