

Product Datasheet

DRD2 Antibody (orb1816924)

Description	DRD2 Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human DRD2
Target	DRD2
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	51 kDa
Uniprot ID	P14416
Expiration Date	12 months from date of receipt.