

Product Datasheet

STAT3 Antibody (orb1817003)

Description	STAT3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human STAT3 around the phosphorylation site of Ser727 PM(p-S)PR
Target	STAT3
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	88 kDa
Uniprot ID	P40763
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.