

Product Datasheet

Recombinant human IL-8 protein, His (HEK293) (orb1817209)

Description	Recombinant human IL-8 protein, His (HEK293)
Endotoxins	Not analyzed
Target	CXCL8
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concent
Protein Sequence	28-99/99
Purity	> 90% as determined by SDS- PAGE
MW	9.9 kDa
Uniprot ID	P10145
Expiration Date	6 months from date of receipt.