

Product Datasheet

Human ORAI2(ORAI Calcium Release-Activated Calcium Modulator 2) ELISA Kit (orb1817345)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.32-20 ng/mL
Note	For research use only
Sample Types	serum, plasma, urine, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	ORAI2
Sensitivity	0.13 ng/mL
Expiration Date	Please enquire.

