

## Product Datasheet

### EasyStep Human OC(Osteocalcin) ELISA Kit (orb1817400)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>      | Human                                                                                                   |
| <b>Range</b>           | 0.23-15 ng/mL                                                                                           |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Sample Types</b>    | Serum                                                                                                   |
| <b>Sensitivity</b>     | 0.14 ng/mL                                                                                              |
| <b>Expiration Date</b> | Please enquire.                                                                                         |

