

Product Datasheet

Anti-TrkB / NTRK2 Reference Antibody (Rinat patent anti-TrkB) (orb1817422)

Description	Anti-TrkB / NTRK2 Reference Antibody (Rinat patent anti-TrkB)
Reactivity	Human
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	TrkB / NTRK2
Isotype	IgG2SA
Clonality	Monoclonal
MW	146.38
Uniprot ID	Q16620
Expiration Date	12 months from date of receipt.