

Product Datasheet

GABA A Receptor Antibody (FITC) (orb181747)

Description

Mouse Anti-Mouse GABA-A Receptor Beta

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	AM, ICC, IHC, WB
Immunogen	Fusion protein amino acids 370-433 of mouse GABA-A-R-Beta 3
Target	GABA A Receptor
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-339 was sufficient for detection of Beta 3 GABA receptor in 10 µg of rat brain lysate by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	55kDa
Uniprot ID	P63080
NCBI	NP_032097.1
Entrez	14402
Dilution Range	WB (1:1000), IHC-P (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.