

Product Datasheet

Anti-DCSTAMP Reference Antibody (U.Rochester patent anti-DC-STAMP) (orb1817608)

Description	Anti-DCSTAMP Reference Antibody (U.Rochester patent anti-DC-STAMP)
Reactivity	Human
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	DCSTAMP
Isotype	IgG1
Clonality	Monoclonal
MW	146.2
Uniprot ID	Q9H295
Expiration Date	12 months from date of receipt.