

Product Datasheet

Anti-CTLA-4 / CD152 Reference Antibody (ipilimumab) (orb1817628)

Description	Anti-CTLA-4 / CD152 Reference Antibody (ipilimumab)
Reactivity	Human, Monkey, Mouse
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	CTLA-4 / CD152
Isotype	IgG1
Clonality	Monoclonal
MW	145.38
Uniprot ID	P16410
Expiration Date	12 months from date of receipt.