



Product Datasheet

EGFP antibody (orb181782)


Description	EGFP antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant EGFP protein (orb84840) was used to make this antibody
Target	EGFP
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	27 kDa
Dilution Range	WB: 500-1500
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant EGF...



WB analysis of recombinant EGF antibody...



Western blot analysis of recombinant EGF...