

Product Datasheet

AFP/Alpha 1 Fetoprotein Monoclonal Antibody (orb18179)

Description

AFP/Alpha 1 Fetoprotein Monoclonal Antibody



Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC

Applications

Immunohistochemical analysis of paraffin...

Immunogen

Human alpha-fetoprotein.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW

68386 MW

Uniprot ID

P02773

Sensitivity

> 5000 cells

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 2-4µg/ml, Human, By Heat

Expiration Date

12 months from date of receipt.