



Product Datasheet

FOXP3 antibody (FITC) (orb181791)


Description	FOXP3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Target	FOXP3 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Dilution Range	IF/ICC: 1:100
Expiration Date	12 months from date of receipt.



IF image of mouse skin tissue using anti...



IF analysis of rat skin tissue using FOX...



Immunofluorescence analysis of rat skin

...