

Product Datasheet

Anti-TNFRSF10A / TRAILR1 / CD261 Reference Antibody (mapatumumab) (orb1817952)

Description	Anti-TNFRSF10A / TRAILR1 / CD261 Reference Antibody (mapatumumab)
Reactivity	Human
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	TNFRSF10A / TRAILR1 / CD261
Isotype	IgG2SA
Clonality	Monoclonal
MW	145 Da
Uniprot ID	000220
Expiration Date	12 months from date of receipt.