

Product Datasheet

HSP25/HSP27 antibody (orb181812)

Description

Mouse monoclonal antibody against Hsp27,

Species/Host	Mouse
Reactivity	Bovine, Canine, Guinea pig, Hamster, Human, Mouse, Rat, Sheep
Conjugation	PerCP
Tested Applications	FACS, FC, ICC, IHC, IP, WB
Immunogen	HSP27 peptide
Target	HSP27
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:5000 dilution of SMC-114 was sufficient for detection of HSP27 in 20 µg of HeLa cell lysate by ECL immunoblot analysis.
Isotype	IgG1 Kappa
Clonality	Monoclonal
MW	27kDa
Uniprot ID	P04792
NCBI	NP_001532.1
Entrez	3315
Dilution Range	WB (1:5000), ICC/IF (1:200), IHC (1:100)
Expiration Date	12 months from date of receipt.