

Product Datasheet

Anti-PDCD1 / PD-1 / CD279 Reference Antibody (camrelizumab) (orb1818318)

Description	Anti-PDCD1 / PD-1 / CD279 Reference Antibody (camrelizumab)
Reactivity	Human, Monkey
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	PDCD1 / PD-1 / CD279
Isotype	IgG4SP
Clonality	Monoclonal
MW	145.76
Uniprot ID	Q15116
Expiration Date	12 months from date of receipt.