

## Product Datasheet

### Anti-IL-3Ra / CD123 Reference Antibody (talacotuzumab) (orb1818422)

|                            |                                                                                                                                  |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Anti-IL-3Ra / CD123 Reference Antibody (talacotuzumab)                                                                           |
| <b>Reactivity</b>          | Human, Monkey                                                                                                                    |
| <b>Endotoxins</b>          | < 0.001EU/μg,determined by LAL method.                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                     |
| <b>Tested Applications</b> | ELISA, FA, FACS, Kinetics                                                                                                        |
| <b>Form/Appearance</b>     | Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column. |
| <b>Concentration</b>       | 1mg/ml                                                                                                                           |
| <b>Note</b>                | For research use only                                                                                                            |
| <b>Protein Sequence</b>    | IL-3Ra / CD123                                                                                                                   |
| <b>Isotype</b>             | IgG1                                                                                                                             |
| <b>Clonality</b>           | Monoclonal                                                                                                                       |
| <b>MW</b>                  | 147.42                                                                                                                           |
| <b>Uniprot ID</b>          | <b>P26951</b>                                                                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                  |