

Product Datasheet

Anti-ERBB1 / EGFR / HER1 Reference Antibody (necitumumab) (orb1818454)

Description	Anti-ERBB1 / EGFR / HER1 Reference Antibody (necitumumab)
Reactivity	Human, Monkey, Mouse
Endotoxins	< 0.001EU/μg,determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Form/Appearance	Purified monoclonal antibody supplied in PBS, pH 6.0, without preservative.This antibody is purified through a protein A column.
Concentration	1mg/ml
Note	For research use only
Protein Sequence	ERBB1 / EGFR / HER1
Isotype	IgG1
Clonality	Monoclonal
MW	144.78
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.