



Product Datasheet

Calbindin-D Calb1 Antibody (Monoclonal, CB-955) (orb18185)


Description

Calbindin-D Calb1 Antibody (Monoclonal, CB-


Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Bovine kidney calbindin-D
Form/Appearance	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	29994 MW
Uniprot ID	P07171
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 1-2µg/ml, Human, bovine, mouse, rat, By Heat Western blot, 0.5-1µg/ml, Human, bovine, mouse, ratImmunofluorescence, 2µg/ml, Mouse
Expiration Date	12 months from date of receipt.



Western Blot of rat tissue lysate Calbin...



IHC-P of rat brain tissue (Calbindin-D a...



Immunohistochemical analysis of paraffin...