

Product Datasheet

Calcineurin Alpha Ppp3r1 Antibody (Monoclonal, CN-A1) (orb18186)

Description

Calcineurin Alpha Ppp3r1 Antibody (Monoclonal,

Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Bovine brain calcineurin.
Form/Appearance	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 ug/ml.
Storage	At -20°C for one year and at 4°C for one month from date of receipt. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Other applications have not been tested. Optimal dilutions should be determined by end users. . Adding 1 ml of PBS buffer will yield a concentration of 100 ug/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	19300 MW
Uniprot ID	P63100
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, bovine, rat, By Heat Western blot, 0.25-0.5µg/ml, Human, bovine, rat
Expiration Date	12 months from date of receipt.

