

Product Datasheet

Trap1 antibody (orb181865)

Description	Mouse monoclonal antibody against Trap1, conjugated to Biotin.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, WB
Immunogen	Purified recombinant TRAP1
Target	TRAP1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-207 was sufficient for detection of Trap-1/HSP75 in 20 µg of Human A431 lysate by ECL immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Isotype	IgG1 Kappa
Clonality	Monoclonal
MW	75kDa
Uniprot ID	Q12931
NCBI	NP_057376.2
Entrez	10131
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.