


Product Datasheet

Alginate Antibody (orb181890)

Description	Mouse Anti- Alginate Monoclonal IgG1 Kappa
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Sodium Alginate conjugated to KLH
Target	Alginate
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-209 was sufficient for detection of 150 ng of alginate-conjugated BSA by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1 Kappa
Clonality	Monoclonal
Dilution Range	WB (1:500)
Expiration Date	12 months from date of receipt.



Western Blot analysis of all bsa-alginat...