

Product Datasheet

**Human IGFBP-4 Protein, Fc Tag (MALS verified)
(orb1818947)**

Description

Insulin-like growth factor-binding protein 4

which...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IGFBP-4

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;
-70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 52.4 kDa. The protein migrates as 60-65 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source


Human IGFBP-4 Protein, Fc Tag (orb1818947) is expressed from human 293 cells (HEK293). It contains AA Asp 22 - Glu 258 (Accession # P22692-1).

MW

52.4 kDa

Expiration Date

6 months from date of receipt.



Immobilized
Human IGF-
I, His Tag,
premiu...