

Product Datasheet

Human Activin RIB & Activin RIIB Protein, Fc Tag&Fc Tag (orb1818950)

Description

Human Activin RIB & Activin RIIB

H5253) is expressed from human 293 cell...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Activin RIB & Activin RIIB

Preservatives

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH 7.5

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc | C-Fc

Note

For research use only

Application notes

This protein . The protein has a calculated MW of 37 kDa & 39.1 kDa. The protein migrates as 80-95 kDa and > 150 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human Activin RIB & Activin RIIB Protein, Fc Tag&Fc Tag (orb1818950) is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Glu 126 & Gly 20 - Pro 133 (Accession # P36896-1 & Q13705-1).