

Product Datasheet

**Influenza A [A/Thailand/1(KAN-1)/2004(H5N1)]
Neuraminidase (NA) Protein, His Tag [Out of
stock] (orb1818954)**

Description

Influenza A [A/Thailand/1(KAN-

Protein, His Tag (HA1-V5245) is
expr...

Reactivity

Virus

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Neuraminidase (NA)

Preservatives

PBS, pH 7.4, 300 mM NaCl

Form/Appearance

Liquid

Storage

The product MUST be stored at -
70°C or lower upon receipt;. -70°C
for 3 months under sterile
conditions

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine
tag at the N-terminus. The protein
has a calculated MW of 53.2 kDa.
The protein migrates as 65 kDa
when calibrated against Star Ribbon
Pre-stained Protein Marker under
reducing (R) condition (SDS-PAGE)
due to glycosylation.

Purity

90%

Source

Influenza A [A/Thailand/1(KAN-
1)/2004(H5N1)] Neuraminidase (NA)
Protein, His Tag (orb1818954) is
expressed from human 293 cells
(HEK293). It contains AA His 36 -
Lys 449 (Accession #
H8PF47_9INFA, GISAID).

MW

53.2 kDa

Uniprot ID
H8PF47
Expiration Date

6 months from date of receipt.