

Product Datasheet

Human CXCL14 Protein, His Tag (orb1818955)

Description

Human CXCL14 Protein, His Tag is

	AA Ser 35 - Glu 111 (Acce...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CXCL14
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 11.4 kDa. The protein migrates as 15 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).
Purity	90%
Source	Human CXCL14 Protein, His Tag (orb1818955) is expressed from E. coli cells. It contains AA Ser 35 - Glu 111 (Accession # O95715).
MW	11.4 kDa
Uniprot ID	O95715
Expiration Date	6 months from date of receipt.