

Product Datasheet

Mouse IL-28B Recombinant Protein (orb1820018)

Description	The Mouse IL-28B yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-28B
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	DPVPRATRLP VEAKDCHIAQ FKSLSPKELQ AFKKAKGAIE KRLLEKDMRC SSSLISRAWD LKQLQVQERP KALQAEVALT LKVWENINDS ALTTILGQPL HTLSHIHSQL QTCTQLQATA EPKPPSRRLS RWLHRLQEAQ SKETPGCLED SVTSNLFQLL LRDLKCVASG DQCV (174)
Purity	98%
Source	Yeast
MW	19.6 kDa
Entrez	338374
Expiration Date	6 months from date of receipt.