

## Product Datasheet

### Mouse IL-28B Recombinant Protein (orb1820020)

|                         |                                                                                                                                                                                                                       |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | The Mouse IL-28B yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...                                                                                                               |
| <b>Target</b>           | IL-28B                                                                                                                                                                                                                |
| <b>Form/Appearance</b>  | Lyophilized                                                                                                                                                                                                           |
| <b>Storage</b>          | -20°C                                                                                                                                                                                                                 |
| <b>Note</b>             | For research use only                                                                                                                                                                                                 |
| <b>Protein Sequence</b> | DPVPRATRLP VEAKDCHIAQ FKSLSPELQ<br>AFKKAKGAIE KRLLEKDMRC SSSLISRAWD<br>LKQLQVQERP KALQAEVALT<br>LKVWENINDS ALTTILGQPL HTLSHIHSQL<br>QTCTQLQATA EPKPPSRRLS<br>RWLHRLQEAQ SKETPGCLED<br>SVTSNLFQLL LRDLCVASG DQCV (174) |
| <b>Purity</b>           | 98%                                                                                                                                                                                                                   |
| <b>Source</b>           | Yeast                                                                                                                                                                                                                 |
| <b>MW</b>               | 19.6 kDa                                                                                                                                                                                                              |
| <b>Entrez</b>           | <b>338374</b>                                                                                                                                                                                                         |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                                                                                                                                        |