

Product Datasheet

Rabbit IL-5 Recombinant Protein (orb1820147)

Description	The Rabbit IL-5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MATEIRMSTV VKETLTLLST YQSLLIGNET LMIPVPVHKN HHLCEETFR GVDTLKAQIV QGEAMDNLQ NLYLIKKYID LQKKKCGEER RGVKHFLDYL QEFLGVINTE WTMES (115)
Purity	98%
Source	Yeast
MW	13.4 kDa
Entrez	100358075
Expiration Date	6 months from date of receipt.