

Product Datasheet

GABA-A Receptor Beta 3 antibody (orb182044)

Description

Mouse monoclonal antibody against

Biotin.

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Biotin

Tested

AM, ICC, IHC, WB

Applications
Immunogen

Fusion protein amino acids 370-433 of mouse GABA-A-R-Beta 3

Target

GABA A Receptor

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-339 was sufficient for detection of Beta 3 GABA receptor in 10 µg of rat brain lysate by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

55kDa

Uniprot ID

P63080

NCBI

NP_032097.1

Entrez

14402

Dilution Range

WB (1:1000), IHC (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.