

Product Datasheet

Guinea Pig TNF alpha Recombinant Protein (orb1820758)

Description	The Guinea Pig TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SASQNDNDKP VAHVVANQQA EEELQWLSKR ANALLANGMG LSDNQLVVPS DGLYLIYSQV LFKGQGCPSY LLLTHTVSRL AVSYPEKVNL LSAIKSPCQK ETPEGAERKP WYEPIYLGGV FQLQKGDRLS AEVNLPQYLD FADSGQIYFG VIAL (154)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	100135630
Expiration Date	6 months from date of receipt.