

Product Datasheet

Guinea Pig IL-8 Recombinant Protein (orb1820759)

Description	The Guinea Pig IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MVVTKLVSEL RCQCIKIHTT PFHPKFIKEL KVIESGPRCA NSEIIVKLSD NRQLCLDPKK KVVQDVVSMF LKRTESQDS (79)
Purity	98%
Source	Yeast
MW	9.1 kDa
Entrez	100379599
Expiration Date	6 months from date of receipt.