


Product Datasheet

Collagen, Type IV COL4A1 Antibody (Monoclonal, COL-94) (orb18208)

Description

Collagen, Type IV COL4A1 Antibody


Immunohistochemical
analysis of paraffin...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC

Applications
Immunogen

Human collagen type IV.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a
concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of
receipt. After reconstitution, at 4°C for one
month. It can also be aliquotted and stored
frozen at -20°C for six months. Avoid
repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested.
Optimal dilutions should be determined by
end users. . Add 1ml of PBS buffer will yield
a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

160615 MW

Uniprot ID
P02462
Dilution Range

Immunohistochemistry (Paraffin-embedded
Section), 0.5-1µg/ml, Human, By Heat
Immunohistochemistry (Frozen Section),
0.5-1µg/ml, Human

Expiration Date

12 months from date of receipt.