

Product Datasheet

Canine IL-6 Recombinant Protein (orb1820872)

Description	The Canine IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPTPGPLAGDSKDDATSNLSLPLTSANKVEELIKYILGKISALRKEMCDKFNKCEDS KEALAENNLHLPKLEGKDGCFQSGFNQETCLTRITTGLVEFQLHLNILQNNYEGD KENVKSVHMSTKILVQMLKSKVKNQDEVTTTPDPTDASLQAILQSQDEWLKHTTI HLILRSLEDFLQFSLRAVRIM
Purity	98%
Source	Yeast
MW	21.0 kDa
Entrez	403985
Expiration Date	6 months from date of receipt.