

Product Datasheet

GABA-B Receptor 2 antibody (orb182095)

Description

Mouse monoclonal antibody against GABA-

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC, WB
Immunogen	Fusion protein amino acids 861-912 of rat GABA(B)R2
Target	GABA B Receptor 2
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-402 was sufficient for detection of GABA(B)R2 in 20 µg of rat brain membrane lysate and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	105kDa
Uniprot ID	O88871
NCBI	NP_113990.1
Entrez	83633
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.