

Product Datasheet

Bovine IL-8 Recombinant Protein (orb1820953)

Description	The Bovine IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVLSRMSTEL RCQCIKTHST PFHPKFIKEL RVIESGPHCE NSEIIVKLTN GNEVCLNPKE KVVQKVVQVF VKRAEKQDP
Purity	98%
Source	Yeast
MW	9.1 kDa
Entrez	280828
Expiration Date	6 months from date of receipt.