



Product Datasheet

CD79a Antibody (FITC) (orb1821019)


Description	CD79a (Ig alpha, MB1) forms disulfide-linked heterodimer with CD79b (Ig beta). They both are transme...
Reactivity	Bovine, Canine, Equine, Gallus, Guinea pig, Human, Mouse, Porcine, Primate, Rabbit, Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Synthetic peptide corresponding to C terminal amino acids 208-222 of human CD79a
Target	CD79a
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests. Intracellular staining.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	P11912
Entrez	973
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
intracellular
...



Separation
of human
CD19
positive
CD79a ...



Flow
cytometry
intracellular
staining pa...