



Product Datasheet

CD151 Antibody (APC) (orb1821028)

Description	CD151, also known as PETA-3 (platelet-endothelial tetraspan antigen), is a four-pass transmembrane g...
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Human epidermoid carcinoma cell line Hep-3
Target	CD151
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	P48509
Entrez	977
Expiration Date	12 months from date of receipt.



Separation of human thrombocytes (red-fi...



Flow cytometry surface staining pattern ...