

Product Datasheet

CHRM3 / M3 Antibody (C-Terminus) (orb1821030)

Description	M3 antibody is an unconjugated rabbit polyclonal antibody to M3 (CHRM3) (C-Terminus) from human. It ...
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 19 amino acid peptide from C-Terminus of human CHRM3 / M3. Percent identity with other species by BLAST analysis: Human, Gorilla, Orangutan, Gibbon, Monkey (100%); Chimpanzee, Mouse, Rat, Sheep, Hamster, Elephant, Dog, Bovine, Horse, Rabbit, Pig, Guinea pig (95%); Opossum, Turkey, Chicken, Lizard (89%).
Preservatives	PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 - 30 µg/ml)
Expiration Date	12 months from date of receipt.