

Product Datasheet

**Lamin A+C Monoclonal Antibody [JOL2] -
Fluorescein (FITC) / 50ul (orb1821243)**

Description

Lamin A/C Monoclonal Antibody:

Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	Recombinant human lamin A
Target	This antibody reacts with both recombinant and native Lamin A and C in humans
Preservatives	Purified antibody conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	p02545
Entrez	103223838, 4000
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at: WB 1:1000, IF 1:100, Flow cytometry 1:1000
Expiration Date	12 months from date of receipt.