

Product Datasheet

Cytokeratin 18 Antibody (orb1821350)

Description	Cytokeratin 18 Monoclonal Antibody [C-04]
Species/Host	Mouse
Reactivity	Human, Mammal
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	Tissue, cells or virus corresponding to Cytokeratin 18. Cytoskeleton preparation of epidermal carcinoma cell line A431
Target	Cytokeratin 18
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P05783
Entrez	3875
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at WB 0.001µg/ml, ICC/IF 1:50-1:100, Flow Cyt 1:50-1:100
Expiration Date	12 months from date of receipt.